

BandIT is the only scientifically designed device to aid in the relief of tennis elbow, carpal tunnel syndrome and other forms of forearm and wrist pain commonly associated with RSI's without cutting off circulation, without causing discomfort and without inhibiting play or work. **BandIT** also is effective in the workplace and is recommended for computer operators, dentists, reporters, writers, factory workers and others who engage in repetitive manual activities.




Pro Band Sports Industries, Inc.,
PO Box 50509, 1483 East Valley Road Suite 18
Santa Barbara, CA 93150-0509
Tel.: 800-676-1550 FAX : 805-969-3928
E-mail: pbsi@banditusa.com

BandIT Style #89 Made in U.S.A.

U.S. Patents 6077241, 5441058, 5642739, 5865782, 3750960
Other U.S. Patents & Foreign Patents Pending

©2013 Pro Band Sports Industries

Instructions for Use



Designed by World Renowned
Orthopaedic Surgeon and
Sports Medicine Specialist,
Donald Fared, M.D.



Left arm placement



Right arm placement



Tightening the strap



Place BandIT on the forearm two finger widths below the elbow crease.

Place **BandIT**[®] on the forearm as shown above. Tighten the strap, pulling it from the underside of your arm towards you. Tighten the strap no more than is necessary to keep **BandIT** in place during the activity that aggravates your condition (i.e., make it tighter for tennis, looser for typing). When worn, the strap between the two plates should be equal in length on each side. The plate with the logo sits on top of the arm; the other, under the arm. To achieve this, pull the bottom plate (the plate without **BandIT**'s logo) along the strap. The strap is *not* glued to the bottom plate. For a more customized fit, the plastic plates are bendable and

may be gently *squeezed* for a smaller arm or *flattened* for a larger arm.

Washing: *DO NOT MACHINE WASH OR DRY.* Hand wash with water and mild soap. Compress pads and rinse thoroughly. With a towel, pat dry any excess water.

Caution: For persistent pain, immediately contact your health care professional. If excessive moisture due to physical exertion builds under the **BandIT**, either place a loose fabric sleeve under it or remove it after activity to allow skin and **BandIT** to dry. *DO NOT USE OVER OILS OR LOTIONS.*